

Rabbit Anti-Cat IgG Heavy & Light Chain antibody {Alexa Fluor 555} (STJS000392)

STJS000392

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 555-conjugated rabbit polyclonal anti-Cat IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA, IF, ICC, FC
Host/Source	Rabbit
Reactivity	Cat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Alexa Fluor 555
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000
	IF 1:200-1:1000
	FCM 1:2000-1:4000
	ELISA Use at an assay dependent concentration
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence