

Rabbit Anti-Human IgG+IgM+IgA antibody {Alexa Fluor 555} (STJS000403)

STJS000403

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 555-conjugated rabbit polyclonal anti-Human IgG+IgM+IgA secondary antibody. For use in most research applications.
Applications	ELISA, IF, ICC, FC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Alexa Fluor 555
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
Isotype	IgG+IgM+IgA
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	