

Rabbit Anti-Monkey IgG Heavy & Light Chain antibody {Alexa Fluor 555} (STJS000405)

STJS000405

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | Alexa Fluor 555-conjugated rabbit polyclonal anti-Monkey IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| Applications | ELISA, IF, ICC, FC |
| Host/Source | Rabbit |
| Reactivity | Monkey |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Alexa Fluor 555 |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range | IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration |
| Formulation | 0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |