

## Mouse Anti-Monkey IgG Heavy & Light Chain antibody {Alexa Fluor 555} (STJS000417)

STJS000417

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Alexa Fluor 555-conjugated mouse monoclonal anti-Monkey IgG Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA, IF, ICC, FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Alexa Fluor 555
<b>Purification</b>	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration
<b>Formulation</b>	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**