

Goat Anti-Mouse IgM Mu Chain antibody {Alexa Fluor 568} (STJS000431)

STJS000431

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 568-conjugated goat polyclonal anti-Mouse IgM Mu Chain secondary antibody. For use in most research applications.
Applications	ELISA, IF, ICC, FC
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Alexa Fluor 568
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	