

## Rabbit Anti-Human IgG+IgM+IgA antibody {Alexa Fluor 680} (STJS000661)

STJS000661

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Alexa Fluor 680-conjugated rabbit polyclonal anti-Human IgG+IgM+IgA secondary antibody. For use in most research applications.
<b>Applications</b>	IF, ICC, FC, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Alexa Fluor 680
<b>Purification</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	Fluorescence 1:50-1:1000 LICOR 1:5000-20000
<b>Formulation</b>	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
<b>Isotype</b>	IgG+IgM+IgA
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	