

Goat Anti-Human IgG Gamma chain antibody {Alexa Fluor 790} (STJS000728)

STJS000728

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Secondary antibodies |
| Short Description | Alexa Fluor 790-conjugated goat polyclonal anti-Human IgG Gamma chain secondary antibody. For use in most research applications. |
| Applications | ELISA, IF, ICC, FC |
| Host/Source | Goat |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Alexa Fluor 790 |
| Purification | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Dilution Range | IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration |
| Formulation | 0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |