

## Goat Anti-Mouse IgG Heavy & Light Chain antibody {Alexa Fluor 790} (STJS000731)

STJS000731

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Alexa Fluor 790-conjugated goat polyclonal anti-Mouse IgG Heavy & Light Chain secondary antibody. For use in most research applications.
<b>Applications</b>	ELISA, IF, ICC, FC
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Alexa Fluor 790
<b>Purification</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Dilution Range</b>	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration
<b>Formulation</b>	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	