

Rabbit Anti-Donkey IgG Heavy & Light Chain antibody {Alexa Fluor 790} (STJS000739)

STJS000739

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 790-conjugated rabbit polyclonal anti-Donkey IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	ELISA, IF, ICC, FC
Host/Source	Rabbit
Reactivity	Donkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Alexa Fluor 790
Purification	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	IHC 1:200-1:1000 IF 1:200-1:1000 FCM 1:2000-1:4000 ELISA Use at an assay dependent concentration
Formulation	0.01M PBS, pH 7.2, 1% BSA, 50% Glycerol, 0.02% Sodium Azide
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	