

Rabbit Anti-Goat IgG antibody {DyLight 405} (STJS000771)

STJS000771

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	DyLight 405-conjugated rabbit polyclonal anti-Goat IgG secondary antibody. For use in most research applications.
Applications	IF, ICC, FC
Host/Source	Rabbit
Reactivity	Goat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	DyLight 405
Purification	Affinity purified using solid phase Goat IgG (H&L) with finally > 95% purity based on SDS-PAGE
Dilution Range	
Formulation	PBS, pH 7.4, 0.02% Sodium Azide, 1% BSA as stabilizer and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	