

Goat Anti-Rabbit IgG antibody {Dylight 594} (STJS000781)

STJS000781

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Dylight 594-conjugated goat polyclonal anti-Rabbit IgG secondary antibody. For use in most research applications.
Applications	IF, ICC, FC
Host/Source	Goat
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Dylight 594
Purification	Affinity purified using solid phase goat IgG (H&L) with finally > 95% purity based on SDS-PAGE
Dilution Range	
Formulation	PBS, pH 7.4, 0.02% Sodium Azide, 1% BSA as stabilizer and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence