

## Goat Anti-Mouse IgG antibody {Dylight 800} (STJS000790)

STJS000790

### GENERAL INFORMATION

<b>Product Type</b>	Secondary antibodies
<b>Short Description</b>	Dylight 800-conjugated goat polyclonal anti-Mouse IgG secondary antibody. For use in most research applications.
<b>Applications</b>	IF, ICC, FC, WB
<b>Host/Source</b>	Goat
<b>Reactivity</b>	Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Dylight 800
<b>Purification</b>	Affinity purified using solid phase Mouse IgG (H&L) with finally > 95% purity based on SDS-PAGE
<b>Dilution Range</b>	Fluorescence 1:50-1:1000 LICOR 1:5000-20000
<b>Formulation</b>	PBS, pH 7.4, 0.02% Sodium Azide, 1% BSA as stabilizer and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	