

Goat Anti-Mouse Mouse antibody {HRP} (STJS000792)

STJS000792

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated goat polyclonal anti-Mouse Mouse secondary antibody. For use in most research applications.
Applications	ELISA, WB, IHC
Host/Source	Goat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	HRP
Purification	IgG fraction purified with Protein A affinity chromatography
Dilution Range	ELISA: 0.1 µg/ml, WB: 0.1 to 0.2 µg/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 50% glycerole
Isotype	IgG
Storage	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	purified mouse IgG fraction. mouse IgG fraction
Immunogen Region	
Specificity	Reacts with mouse IgG-s
Immunogen Sequence	