

Rabbit Anti-Chicken Chicken antibody {HRP} (STJS000795)

STJS000795

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	HRP-conjugated rabbit polyclonal anti-Chicken Chicken secondary antibody. For use in most research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Chicken

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	HRP
Purification	IgG fraction purified by IgY affinity chromatography
Dilution Range	ELISA: 0.03 to 0.05 µg/ml, WB: 0.1 to 0.2 µg/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 50% glycerol
Isotype	IgG
Storage	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified IgY fraction. chicken IgY
Immunogen Region	
Specificity	Reacts with chicken IgY antibodies
Immunogen Sequence	