

Rabbit Anti-Mouse IgG (H&L) antibody {Alexa Fluor 555} (STJS000806)

STJS000806

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 555-conjugated rabbit polyclonal anti-Mouse IgG H&L secondary antibody. For use in most research applications.
Applications	
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Alexa Fluor 555
Purification	Affinity purified
Dilution Range	Fluorescence microscopy: 1-2 ug/mL, Flow cytometry: 1 ug/10 ⁶ cells, Near-infrared western detection: 50-100 ng/mL
Formulation	Pink solution in pH~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.
Isotype	IgG
Storage	Store at 4°C for several months. Protect from light. For longer storage, divide the conjugate into small aliquots and freeze at -20°C.
Instruction	Avoid repeated freezing and thawing.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	H&L
Immunogen	H&L
Region	
Specificity	
Immunogen	
Sequence	