

Donkey Anti-Donkey IgG Highly Cross-Adsorbed (H&L) antibody {Alexa Fluor 647} (STJS000846)

STJS000846

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Alexa Fluor 647-conjugated donkey polyclonal anti-Donkey IgG Highly Cross-Adsorbed H&L secondary antibody. For use in most research applications.
Applications	
Host/Source	Donkey
Reactivity	Donkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Alexa Fluor 647
Purification	Affinity purified
Dilution Range	Fluorescence microscopy: 1-2 ug/mL. Flow cytometry: 1 ug/10 ⁶ cells. Near-infrared western detection: 50-100 ng/mL
Formulation	Blue solution in pH~7.4 PBS containing 50% glycerol, 2 mg/ml bovine serum albumin (IgG-free and protease-free) and 0.05% sodium azide.
Isotype	IgG
Storage	Store at 4°C for several months. Protect from light. For longer storage, divide the conjugate into small aliquots and freeze at -20°C.
Instruction	Avoid repeated freezing and thawing.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	H&L
Immunogen	H&L
Region	
Specificity	
Immunogen	
Sequence	