

Anti-SCG3 antibody (199-468) (STJ110109)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 53kDa Observed Mw: 53kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29106
Gene Symbol	SCG3
UniProt ID	SCG3_HUMAN
Immunogen Region	199-468
Immunogen Sequence	KLISKEANNYEEDPNKPTSW TENQAGKIPEKVTPMAAIQD GLAKGENDETVSNTLTLTNG LERRTKTYSEDNFEELQYFP NFYALLKSIDSEKEAKEKET LITIMKTLIDFVKMMVKYGT ISPEEGVSYLENLDEMIALQ TKNKLEKNATDNISKLFAP SEKSHEETDSTKEEAAKMEK EYGSLKDKSTKDDNSNPGGKT DEPKGKTEAYLEAIRKNIEW LKKHDKKGNKEDYDLKMR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 199-468 of human SCG3 (NP_037375.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using SCG3 antibody (STJ110109) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.