

Anti-ZNF350 antibody (285-532) (STJ110139)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 60kDa Observed Mw: 60kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59348
Gene Symbol	ZNF350
UniProt ID	ZNF350_HUMAN
Immunogen Region	285-532
Immunogen Sequence	TGEKPYICSECGKGFQKGN LIVHQRHTGEKPYICNECG KGFIQKTCLIAHQRFHTGKT PFVCSECGKSCSQKSLIKH QRIHTGEKPFECSECGKAFS TKQKLIVHQRHTHTGERPYGC NECGKAFAYMSCLVKHKRIH TREKQEAARKVENPPAERHSS LHTSDVMQEKNSANGATTQV PSVAPQTSLSNISGLLANRNV VLVGQPVVRCASGDNRGFA QDRNLVNAVNVVPSVINY
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 285-532 of human ZNF350 (NP_067645.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of Jurkat, using ZNF350 Rabbit polyclonal antibody (STJ110139) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.