

Anti-ORC6 antibody (1-252) (STJ110273)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 28kDa Observed Mw: 28kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23594
Gene Symbol	ORC6
UniProt ID	ORC6_HUMAN
Immunogen Region	1-252
Immunogen Sequence	MGSELIGRLAPRLGLAEPDM LRKAE EYLR LSRVKCVGLSA RTTETSSAVMCLDLAASWMK CPLDRAYLIKLSGLNKETYQ SCLKSFEC LLGLNSNIGIRD LAVQFSCIEAVNMASKILKS YESSLPQTQQVDLDLSRPLF TSAALLSACKILKLKVDKNK MVATSGVKAIFDR LCKQLE KIGQQVDREP GDVATPPRKR KKIVVEAPAKEMEKVEEMPH KPQKDELTQDYEEWKRKI
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-252 of human ORC6 (NP_055136.1).

ADDITIONAL INFORMATION

Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using ORC6 Rabbit pAb (STJ110273) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit Exposure time: 30s.