

Anti-TRMT61B antibody (1-250) (STJ110370)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 53kDa Observed Mw: 52kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55006
Gene Symbol	TRMT61B
UniProt ID	TR61B_HUMAN
Immunogen Region	1-250
Immunogen Sequence	MLMAWCRGPPVLLCLRQGLGT NSFLHGLGQEPFEGARSLCC RSSPRDLRDGEREHEAAQRK APGAESCPSLPLSISDGTG CLSSLENLRLPTLREESSPR ELEDSSGDQGRCPHQQGSE DPSMLSQAQSATEVEERHVS PSCSTSRERPFQAGELILAE TGETKFKKLFRLNNFGLL NSNWGAVPFGKIVGKFPQQI LRSSFGKQYMLRRPALEDYV VLMKRGTAITFPKINMIL
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TRMT61B (NP_060380.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using T
 Secondary antibody: HRP-conjugated Goat anti-
 Rabbit IgG (H+L) (STJS000856) at 1:10000
 dilution.

Lysates/proteins: 25 Mu g per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 5s.