

Anti-ATP5O antibody (24-213) (STJ110402)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 23kDa Observed Mw: 23kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	539
Gene Symbol	ATP5PO
UniProt ID	ATPO_HUMAN
Immunogen Region	24-213
Immunogen Sequence	FAKLVRPPVQVYIEGGRYAT ALYSAASKQNKLEQVEKELL RVAQILKEPKVAASVLNPYV KRSIKVKSLNDITAKERFSP LTTNLLINLLAENGRSLNTQG VVSFAFSTMMSVHRGEVPCTV TSASPLEEATLSELKTVLKS FLSQGGQVLKLEAKTDPSILG GMIVRIGEKEYVDMMSVKTKIQ KLGRAMREIV
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-213 of human ATP5O (NP_001688.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of various lysates using ATP5O Rabbit pAb (STJ110402) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.
 Lysates/proteins: 25 µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit
 Exposure time: 90s.



Immunohistochemistry analysis of paraffin-embedded Rat stomach using ATP5O Rabbit pAb (STJ110402) at dilution of 1:100 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to immunohistochemistry staining.



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer using ATP5O Rabbit pAb (STJ110402) at dilution of 1:100 (40x lens). High pressure antigen retrieval performed with 0.01M Citrate Bufferr (pH 6.0) prior to immunohistochemistry staining.



Immunofluorescence analysis of U2OS cells using ATP5O Rabbit pAb (STJ110402) at dilution of 1:100. Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (STJS001166) at 1:500 dilution. Blue: DAPI for nuclear staining.