

Anti-INADL antibody (450-680) (STJ110774)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 196kDa Observed Mw: 255kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10207
Gene Symbol	PATJ
UniProt ID	INADL_HUMAN
Immunogen Region	450-680
Immunogen Sequence	VRRKTSSTSPLEPPSDRGT VVEPLKPPALFLTGA VETET NV DGEDEEIKERIDTLKNDN IQALEKLEKVPDSPENELKS RWENLLGPDYEVMTLDLTDQ IADDAELQKYSKLLPIHTLR LGVEVDSFDGHHYISSIVSG GPVDTLGLLPEDDELLEVNG MQLYGKSRREAVSFLKEVPP PFTLVCCRRFLDDEASVDEP RRTETSLPETEVDHNMDVNT EEDDDGELALW
Specificity	A synthetic peptide corresponding to a sequence within amino acids 450-680 of human INADL (NP_795352.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using INADL Rabbit pAb (STJ110774) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 10s.