

Anti-GTF3C3 antibody (1-160) (STJ11100359)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 101kDa Observed Mw: 101kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	9330
Gene Symbol	GTF3C3
UniProt ID	TF3C3_HUMAN
Immunogen Region	1-160
Immunogen Sequence	MSGFSPELIDYLEGKISFEE FERRREERKTREKKSLEKEG KLSAEENPDDSEVPSSSSGIN STKSQDKDVNEGETSDGVRK SVHKVFASMLGENEDEDEEEE EEEEEEEEEETPEQPTAGD VFVLEMVLNRETKKMMKEKR PRSKLPRLRGLMGANIRF
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human GTF3C3 (NP_036218.1).

ADDITIONAL INFORMATION


Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---




Immunofluorescence analysis of U2OS cells using GTF3C3 Rabbit polyclonal antibody (STJ11100359) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of H9C2 cells using GTF3C3 Rabbit polyclonal antibody (STJ11100359) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of GTF3C3 in paraffin-embedded Mouse brain using GTF3C3 Rabbit polyclonal antibody (STJ11100359) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of GTF3C3 in paraffin-embedded Rat brain using GTF3C3 Rabbit polyclonal antibody (STJ11100359) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using GTF3C3 Rabbit polyclonal antibody (STJ11100359) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.