

Anti-N6-Methyl-m6A antibody [S4MM] (STJ11100594)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	DB/ELISA
Host / Source	Mouse
Reactivity	Species Independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:500-1:2000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	Chemical compounds corresponding to N6-methyladenosine/m6A.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



The m6A mouse monoclonal antibody (STJ11100594) are tested in Dot Blot against N6-methyladenosine (m6A) and unmodified adenosine (100ng for each oligomer). Oligomer 8- N6-methyladenosine (ATAACTGG-m6A-CCGAATGG) Oligomer 9- unmodified adenosine (ATAACTGGACCGAATGG) Oligomer 10- unmodified adenosine (AAAAAAAAAAAAA-biotin)