

Anti-GTF2F1 antibody (390-517) (STJ11100599)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 58kDa Observed Mw: 74kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2962
Gene Symbol	GTF2F1
UniProt ID	T2FA_HUMAN
Immunogen Region	390-517
Immunogen Sequence	PSAEGGSTSSTLRAAASKLE QGKRVSEMPAAKRLRLDTGP QSLSGKSTPQPPSGKTTTNS GDVQVTEDAVRRYLTRKPM T KDLLKKFQTKKTGLSSEQT VNVLAQILKRLNPERKMIND KMHFSLKE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 390-517 of GTF2F1 (NP_002087.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using GTF2F1 Rabbit polyclonal antibody (STJ11100599) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.