


Anti-CRYZ antibody (1-329) (STJ11100873)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 35kDa Observed Mw: 35kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	1429
	Gene Symbol	CRYZ
	UniProt ID	QOR_HUMAN
	Immunogen Region	1-329
	Immunogen Sequence	MATGQKLMRAVRVFEFGGPE VLKLRSDIAVPIPKDHQVLI KVHACGVNPVETYIRSGTYS RKPLLPTYPGSDVAGVIEAV GDNASAFKKGDRVFTSSTIS GGYAEYALAADHTVYKLPEK LDFKQGAAGIPYFTAYRAL IHSACVKAGESVLVHGASGG VGLAACQIARAYGLKILGTA GTEEGQKIVLQNGAHEVFNH REVNYIDKIKKYVGEKGIDI IIEMLANVNL SKDLSLLSH
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-329 of human CRYZ (NP_001123514.1).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from wild type (WT) and CRYZ knockout (KO) HeLa cells, using [KO Validated] CRYZ Rabbit polyclonal antibody (STJ11100873) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.