

Anti-GABRB2 antibody (300-400) [S9MR] (STJ11101359)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 59kDa Observed Mw: 63kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2561
Gene Symbol	GABRB2
UniProt ID	GABRB2_HUMAN
Immunogen Region	300-400
Immunogen Sequence	PYVKAIDMYLMGCFVFVFMALLEYALVNYIFFGRGPQRQK KAAEKAASANNEKMRLDVNK IFYKDIKQNGTQYRSLWDPTGNLSPTRRRTTNYDFSPLYTMD P
Specificity	A synthetic peptide corresponding to a sequence within amino acids 300-400 of human GABA A Receptor beta 2 (GABRB2) (GABRB2) (P47870).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using GABA A Receptor beta 2 (GABRB2) Rabbit monoclonal antibody (STJ11101359) at 1: 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.