

Anti-DLG3 antibody (464-817) (STJ11101438)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:100-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 90kDa Observed Mw: 110kDa/115kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1741
Gene Symbol	DLG3
UniProt ID	DLG3_HUMAN
Immunogen Region	464-817
Immunogen Sequence	QYRPEEYSRFESKIHDLREQ MMNSSMSSGSGSLRTSEKRS LYVRALFDYDRTRDSCLPSSQ GLSFSYGDILHVINASDDEW WQARLVTPHGESEIQGVIPS KKRVEKKERARLKTVKFHAR TGMIESNRDFFGLSDDYYGA KNLKGQEDAILSYPEVTRQE IHYARPVIIIGPMKDRVND LISEFPHKFGSCVPHTTTPR RDNEVDGQDYHFVVSREQME KDIQDNKFIEAGQFNDNLY
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 464-817 of human DLG3 (NP_066943.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunohistochemistry analysis of DLG3 in paraffin-embedded rat brain using DLG3 Rabbit polyclonal antibody (STJ11101438) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using DLG3 polyclonal antibody (STJ11101438) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.