

Anti-PTPRC antibody (194-390) [S9MR] (STJ11101559)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 IHC-P:1:200-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 147kDa Observed Mw: 180-240kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5788
Gene Symbol	PTPRC
UniProt ID	PTPRC_HUMAN
Immunogen Region	194-390
Immunogen Sequence	SDAYLNASETTTLSPSGSAV ISTTTIATTPSKPTCDEKYA NITVDYLYNKETKLF TAKLN VNENVECGNNTCTNNEVHNL TECKNASVSISHNSCTAPDK TLILDVPPGVEKFQLHDC TQ VEKADTTTICLKWKNIETFTC DTQNITYRFQCGNMIFDNKE IKLENLEPEHEYKCDSEILY NNHKFTNASKIIKTDFG
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 194-390 of human CD45 (NP_002829.3).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human tonsil using CD45 Rabbit monoclonal antibody (STJ11101559) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of Jurkat, using CD45 Rabbit monoclonal antibody (STJ11101559) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.