

## Anti-NFATC1 antibody (844-943) [S1MR] (STJ11101721)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                            |
| Concentration       | Lot specific                                                                                                                                          |
| Conjugation         | Unconjugated                                                                                                                                          |
| Purification        | Affinity purification                                                                                                                                 |
| Dilution Range      | WB:1:500-1:1000<br>ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.                                                                                         |
| Isotype             | IgG                                                                                                                                                   |
| Molecular Weight    | Protein Mw: 101kDa<br>Observed Mw: 110kDa                                                                                                             |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                        |

### TARGET INFORMATION

|                    |                                                                                                             |
|--------------------|-------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">4772</a>                                                                                        |
| Gene Symbol        | <a href="#">NFATC1</a>                                                                                      |
| UniProt ID         | <a href="#">NFAC1_HUMAN</a>                                                                                 |
| Immunogen Region   | 844-943                                                                                                     |
| Immunogen Sequence | PSSPSPPLPPATQEPTCLQP CSPACPPATGRPOHLPSTVR RDESPTAGPRLLPEVHEDGS<br>PNLAPIPVTVKREPEELDQL YLDDVNEIIRNDLSSTSTHS |
| Specificity        | A synthetic peptide corresponding to a sequence within amino acids 844-943 of human NFAT2 (O95644).         |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using NFAT2 antibody (STJ11101721) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.