

Anti-XPO5 antibody (1-200) [S4MR] (STJ11101954)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 136kDa Observed Mw: 136kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	57510
Gene Symbol	XPO5
UniProt ID	XPO5_HUMAN
Immunogen Region	1-200
Immunogen Sequence	MAMDQVNALCEQLVKAVTVM MDPNSTQRYRLEALKFCEEF KEKCPICVPCGLRLAEKTQV AIVRHFGQLQILEHVVKFRWN GMSRLEKVKYLNKNSVMELIAN GTLNILEEENHIKDALSRIV VEMIKREWPHQHPDMLIELD TLSKQGETQTELVMFILLRL AEDVVTFQTLPPQRRRDIQQ TLTQNMERIFSLLNLTQEN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Exportin 5 (XPO5) (Q9HAV4).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human esophageal cancer using Exportin 5 (XPO5) (XPO5) Rabbit monoclonal antibody (STJ11101954) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using Exportin 5 (XPO5) (XPO5) Rabbit monoclonal antibody (STJ11101954) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Western blot analysis of extracts of various cell lines, using Exportin 5 (XPO5) (XPO5) Rabbit monoclonal antibody (STJ11101954) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.