

Anti-RAD50 antibody (1200-1300) [S1MR] (STJ11101961)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 154kDa Observed Mw: 154kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10111
Gene Symbol	RAD50
UniProt ID	RAD50_HUMAN
Immunogen Region	1200-1300
Immunogen Sequence	RCSAGQKVLASLIIRLALAE TFCLNCGIILDEPTTNLDR ENIESLAHALVEIHKRSRQQ RNFQLLVITHDEDFVELLGR SEYVEKFYRIKKNIDQCSEI V
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1200-1300 of human Rad50 (Q92878).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using Rad50 Rabbit monoclonal antibody (STJ11101961) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.