

Anti-RPS6KA1 antibody (1-100) [S3MR] (STJ11102133)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 83kDa Observed Mw: 83kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6195
Gene Symbol	RPS6KA1
UniProt ID	KS6A1_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MPLAQLKEPWPLMELVPLDP ENGQTSGEELQPSKDEGV LKEISITHVVKAGSEKADPS HFELLKVLGQGSFGKVFLVR KVTRPD SGHLYAMKVLKKAT
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human p90Rsk/RSK1/RPS6KA1 (Q15418).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using p90Rsk/RSK1/RPS6KA1 Rabbit monoclonal antibody (STJ11102133) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3min.