

Anti-TSC1 antibody (1065-1164) [S3MR] (STJ11102263)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 130kDa Observed Mw: 150kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7248
Gene Symbol	TSC1
UniProt ID	TSC1_HUMAN
Immunogen Region	1065-1164
Immunogen Sequence	TTMGEASASIPPTTVGSLPSS KSFLGMKARELFRNKSESQC DEDGMTSSLSESLKTELGKD LGVEAKIPLNLDGPHSPPT PDSVQQLHIMDYNETHHEHS
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1065-1164 of human Hamartin (Q92574).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of HT-1080 cells, using TSC1 Rabbit monoclonal antibody (STJ11102263) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 3min.