

Anti-KIF2C antibody (1-100) [S4MR] (STJ11102434)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 81kDa Observed Mw: 85kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11004
Gene Symbol	KIF2C
UniProt ID	KIF2C_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MAMDSSLQARLFPGLAIKIQ RSNGLIHSANVRTVNLEKSC VSVEWAEGGATKGKEIDFDD VAAINPELLQLLPLHPKDNL PLQENVTIQKQKRRSVNSKI
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIF2C/KIF2C/MCAK (Q99661).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using KIF2C/KIF2C/MCAK Rabbit mAb (STJ11102434) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit

Exposure time: 3min.