



Anti-KIFC1 antibody (1-100) [S3MR] (STJ11102673)


GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/IF/ICC/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 74kDa Observed Mw: 74kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3833
	Gene Symbol	KIFC1
	UniProt ID	KIFC1_HUMAN
	Immunogen Region	1-100
	Immunogen Sequence	MDPQRSPLLEVKGNIELKRP LIKAPSQLPLSGSRLKRRPD QMEDGLEPEKKRTRGLGATT KITTSHPRVPSLTTVPQTQG QTTAQKVSKKTGPRCSTAIA
	Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human KIFC1 (Q9BW19).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunofluorescence analysis of U-2 OS cells using KIFC1 Rabbit monoclonal antibody (STJ11102673) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH-3T3 cells using KIFC1 Rabbit monoclonal antibody (STJ11102673) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of paraffin-embedded human appendix using KIFC1 Rabbit monoclonal antibody (STJ11102673) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM Tris/EDTA buffer pH 9.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using KIFC1 antibody (STJ11102673) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.