

Anti-STK36 antibody (250-450) (STJ11102733)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 144kDa Observed Mw: 144kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27148
Gene Symbol	STK36
UniProt ID	STK36_HUMAN
Immunogen Region	250-450
Immunogen Sequence	YHPFIAGHVITIIEPAGPDL GTPFTSRLPPELQVLKDEQA HRLAPKGNQSRILTQAYKRM AEEAMQKKHQTGPALQED KTSKVAPGTAPLPRLGATPQ ESSLLAGILASELKSSWAKS GTGEVPSAPRENRTTPDCER AFPEERPEVLGQRSTDVVDL ENEEPDSDNEWQHLETTEP VPIQLKAPLTLLCNPDFCQR I
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 250-450 of human STK36 (NP_056505.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using STK36 Rabbit polyclonal antibody (STJ11102733) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.