

Anti-NT5E antibody (32-264) (STJ11102804)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 63kDa Observed Mw: 70kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	4907
Gene Symbol	NT5E
UniProt ID	5NTD_HUMAN
Immunogen Region	32-264
Immunogen Sequence	LHTNDVHSRLQTSSEDSSKC VNASRCMGGVARLFTKVQQI RRAEPNVLLLDAGDQYQGTI WFTVYKGA EVAHFMNALRYD AMALGNHEFDNGVEGLIEPL LKEAKFPILSANIKAKGPLA SQISGLYLPYKVLPGDEVV GIVGYTSKETPFLSNPGTNL VFEDITALQPEVDKCLKTLN VNKIIALGHSGFEMDKLIAQ KVRGVDDVVGGHSNTFLYTG NPPSKEVPAGKYK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 32-264 of human CD73/NT5E (NP_002517.1).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of extracts of various cell lines, using CD73/NT5E rabbit polyclonal antibody (STJ11102804) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunohistochemistry analysis of paraffin-embedded rat kidney using CD73/NT5E rabbit polyclonal antibody (STJ11102804) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human placenta using CD73/NT5E rabbit polyclonal antibody (STJ11102804) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse testis using CD73/NT5E rabbit polyclonal antibody (STJ11102804) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.