

Anti-ZNF75A antibody (200-296) (STJ11102827)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 35kDa Observed Mw: 35kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7627
Gene Symbol	ZNF75A
UniProt ID	ZN75A_HUMAN
Immunogen Region	200-296
Immunogen Sequence	WSSDLNKHLLTTHQGIPKYK SWCGKSFSQNTNLHTHQRTH TGEKPFTCHECGKKFSQNSH LIKHRRTHTGEQPYTCSICR RNFSRRSSLLRHQKLHL
Specificity	A synthetic peptide corresponding to a sequence within amino acids 200-296 of human ZNF75A (NP_694573.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from A-549 cells, using ZNF75A Rabbit polyclonal antibody (STJ11102827) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 300s.