

Anti-FAM98B antibody (1-225) (STJ11103367)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 46kDa Observed Mw: 48kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	283742
Gene Symbol	FAM98B
UniProt ID	FA98B_HUMAN
Immunogen Region	1-225
Immunogen Sequence	MRGPEPGPQPTMEGDVLDLTL EALGYKGPLLEEALTKAAE GGLSSPEFSELGIWLGSIQIK SLCNLEESITSAGRDDLESF QLEISGFLKEMACPYSVLIS GDIKDRLLKKEDCLKLLFL STELQASQILQNKHKHNSQL DKNSEVYQEVQAMFDLTLGIP KSTTSDIPHMLNQVESKVKD ILSKVQKNHVGKPLMKMDLN SEQAQLERINDALSCEYEC RRRML
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human FAM98B (NP_775882.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using FAM98B Rabbit polyclonal antibody (STJ11103367) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.