

## Anti-ZNF521 antibody (210-420) (STJ11103408)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 148kDa Observed Mw: 148kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">25925</a>
Gene Symbol	<a href="#">ZNF521</a>
UniProt ID	<a href="#">ZNF521_HUMAN</a>
Immunogen Region	210-420
Immunogen Sequence	GFLSSSSLGHMQVHERNKD GSQSGSRMEDWKMKDTQKCS QCEEGFDFPEDLQKHIAECH PECSPNEDRAALQCVCHEL FVEETSLMNHMEQVHSGEKK NSCSICSESFHTVEELYSHM DSHQQPESCNSNSPSLTVV GYTSVSSSTTPDSNLSVDSST MVEAAPPIPKSRGRKRAAQQ TPDMTGPPSSKQAKVTYSYCIY CNKQLFSSLAIV
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 210-420 of human ZNF521 (NP_056276.1).

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Rat brain, using ZNF521 Rabbit polyclonal antibody (STJ11103408) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25  $\mu$ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.