

Anti-HK1 antibody (20-300) (STJ11103807)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 102kDa Observed Mw: 110kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3098
Gene Symbol	HK1
UniProt ID	H XK1_HUMAN
Immunogen Region	20-300
Immunogen Sequence	KKIDKLYAMRLSDETLIDI MTRFRKEMKNGLSRDFNPTA TVKMLPTFVRSIPDGSEKGD FIALDLGGSSFRILRVQVNH EKNQNVHMESEVYDTPENIV HGSGSQLFDHVAECLGDFME KRKIKDKKLPVGFTFSFPCQ QSKIDEAILITWTKRFKASG VEGADVVKLLNKAIKKRGDY DANIVAVVNDTVGTMTCGY DDQHCEVGLIIGTGTNACYM EELRHIDLVEGDEGRMCIN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human HK1 (NP_000179.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of 293T cells, using HK1 antibody (STJ11103807). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST.