



Anti-ZNF213 antibody (132-459) (STJ11104059)


GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
	Isotype	IgG1
	Molecular Weight	Protein Mw: 51kDa Observed Mw: Refer to figures
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	7760
	Gene Symbol	ZNF213
	UniProt ID	ZN213_HUMAN
	Immunogen Region	132-459
	Immunogen Sequence	RQDVPSEEAPEAAAGRGSQA TGPPTVGARRRPSVPQEQH SHSAQPPALLKEGRPGGETTD TCFVSGVHGPVALGDIPFYF SREEWGTLDPAQRDLFWDIK RENSRTTLGFGGLKGQSEKS LLQEMVPVVGQTGSDVTVS WSPEEAEAWESENRPRAALG PVVGARRGRPPTRRRQFRDL AAEKPHSCGQCGRFRWGS D LARHQRTHTGEKPHKCPECD KSFRSSSDLVRHQGVHTGE
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 132-459 of human ZNF213 (NP_004211.1).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Immunohistochemistry analysis of paraffin-embedded rat testis using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded mouse brain using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human lung cancer using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using ZNF213 Rabbit polyclonal antibody (STJ11104059) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.