

Anti-Um antibody [S4062RM] (STJ11104062)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	DB/ELISA
Host / Source	Rabbit
Reactivity	Species Independent

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	DB:1:500-1:1000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Sequence	Um conjugated to KLH
Specificity	Um

ADDITIONAL INFORMATION

Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Dot-blot analysis of 2-NA-O-methyluridine (Um) and unmodified uridine using 2-NA-O-methyluridine (Um) Rabbit monoclonal antibody (STJ11104062) at 1:1000 dilution. NA 2-NA-O-methyluridine : Biotin-5 NA UGACAACUACAGAC (Um) 3 NA NA NA unmodified uridine : Biotin-5 NA UGACAACUACAGACU3 NA