

Anti-STRADA antibody (25-210) (STJ11104490)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 48kDa Observed Mw: Refer to figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	92335
Gene Symbol	STRADA
UniProt ID	STRAA_HUMAN
Immunogen Region	25-210
Immunogen Sequence	FLPEGGCYELLTVIGKGFED LMTVNLARYKPTGEYVTVRR INLEACSNEMVTFLLQGEHLV SKLFNHPNIVPYRATFIADN ELWVVTSMAYGSAKDLICT HFMDGMNELAIAYILQGVLK ALDYIHMMGYVHRSVKASHI LISVDGKVYLSGLRSNLSMI SHGQRQRVVHDFPKYSVKVL PWLSPE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 25-210 of human STRADA (NP_699166.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from Mouse skeletal muscle, using STRADA Rabbit polyclonal antibody (STJ11104490) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.