

## Anti-JMJD6 antibody (1-170) (STJ11104574)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 46kDa Observed Mw: 46kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">23210</a>
Gene Symbol	<a href="#">JMJD6</a>
UniProt ID	<a href="#">JMJD6_HUMAN</a>
Immunogen Region	1-170
Immunogen Sequence	MNHKSKKRIREAKRSARPEL KDSLDWTRHNYESFSLSPA AVADNVERADALQLSVEEFV ERYERPYPKPVLLNAQEGWS AQEKWTLERLKRKYRNQKFK CGEDNDGYSVKMKMKYYIEY MESTRDDSPLYIFDSSYGEH PKRRKLLLEDYKVPKFFTTDDL FQYAGEKRRP
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human JMJD6 (NP_055982.2).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunohistochemistry analysis of JMJD6 in paraffin-embedded human lung using JMJD6 Rabbit polyclonal antibody (STJ11104574) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of lysates from Mouse liver using JMJD6 Rabbit polyclonal antibody (STJ11104574) at 1:1000 dilution. Secondary antibody:HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25  $\mu$ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection:ECL Basic Kit. Exposure time:90s.