

Anti-NFKB2 antibody (800-900) [S4787RM] (STJ11104787)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 97kDa Observed Mw: 120kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4791
Gene Symbol	NFKB2
UniProt ID	NFKB2_HUMAN
Immunogen Region	800-900
Immunogen Sequence	LRSLVDTYRQTTSPSGSLLR SYELAGGDLA GLLEALSDMG LEEGVRLLRGPETRDKLPST AEVKEDSAYGSQSVEQEAEL LGPPPEPPGGLCHGHPQPQV H
Specificity	A synthetic peptide corresponding to a sequence within amino acids 800-900 of human NF-Kappa B2 (NP_001070962.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates, using NF-Kappa B2 antibody (STJ11104787) at 1:700 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.