

Anti-CD98 antibody (105-526 aa) (STJ11104992)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 58kDa Observed Mw: 75-120kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	17254
Gene Symbol	Slc3a2
UniProt ID	4F2_MOUSE
Immunogen Region	105-526 aa
Immunogen Sequence	ELPVQRWWHKGALYRIGDLQ AFVGRDAGGIAGLKSHLEYL STLKVKGLVLGPIHKNQKDE INETDLKQINPTLGSQEDFK DLLQSAKKKSIHILDTPN YQGQNAWFLPAQADIVATKM KEALSSWLQDGVDFQFRDV GKLMNAPLYLAEWQNTKNL SEDRLLIAGTESSDLQQIVN ILESTSDLLLTSSYLSNSTF TGERTESLVTRFLNATGSQW CSWSVSQAGLLADFIPDHL
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 105-526 of mouse CD98. (NP_032603.3).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates, using CD98 antibody (STJ11104992) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Western blot analysis of various lysates, using CD98 antibody (STJ11104992) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.