

Anti-FCER2 antibody (48-321 aa) [S5017RM] (STJ11105017)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 47kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2208
Gene Symbol	FCER2
UniProt ID	FCER2_HUMAN
Immunogen Region	48-321 aa
Immunogen Sequence	DTTQSLKQLEERAARNVSQV SKNLESHHGDQMAQKSQSTQ ISQELEELRAEQQRLLKSQDL ELSWNLNGLQADLSSFKSQE LNERNEASDLLERLREEVTK LRMEQLQVSSGFVCNTCPEKW INFQRKCYFFGKGTKQWVHA RYACDDMEGQLVSIHSPEEQ DFLTKHASHTGSWIGLRNLD LKGEFIWVDGSHVDYSNAP GEPTSRSQGEDCVMMRGSGR WNDAFCDRKLGAWVCDRLA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 48-321 of human CD23 (NP_001993.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using CD23 antibody (STJ11105017) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.