

Anti-CD7 antibody (26-180 aa) [S5109RM] (STJ11105109)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:5000-1:20000 ELISA: Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 0.05% BSA, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 25kDa Observed Mw: 38-40kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	924
Gene Symbol	CD7
UniProt ID	CD7_HUMAN
Immunogen Region	26-180 aa
Immunogen Sequence	AQEVQQSPHCTTVPVGASVN ITCSTSGGLRGIYLRQLGPQ PQDIIYYEDGVVPTDRRFR GRIDFSGSQDNLITIMHRLQ LSDTGTYTCQAITEVNVYGS GTLVLVTEEQSQGWHRCSDA PPRASALPAPPTGSALPDQP TASALPDPPAASALP
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7. (NP_006128.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using CD7 (STJ11105109) at 1:20000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.